



NovaFount FS-830

One-Step Acid Fountain Solution for Offset Lithography

Description

NovaFount FS-830 is a one-step acid fountain solution for offset lithography.

Directions

Use two to six ounces per gallon as needed. Suggested mix is four ounces per gallon.

Product Need and Usage

This product is formulated for those who need a one-step fountain solution. It is specifically designed for improved wetting and ink/water balance at restart of press for Heidelberg XL, among others.

Features

- * Excellent restart.
- * Excellent ink/water balance.
- * Does not inhibit drying.

Technical Specifications

VOC (lb/gal method 24)	4.52
Specific Gravity	1.050
pH	3.5
Dosage	2 to 6 oz./gal.
Color	Clear, Colorless

Rev. 6/23/2010

Disclaimer: Unless otherwise restricted by applicable law, nothing contained in this literature shall be deemed a representation or warranty of any kind, either expressed or implied. See the Material Safety Data Sheets (MSDS) for this product for safety information prior to use.